



OPERATORY COMPUTER SPECIFICATIONS

	Recommended	BEST Performance
Processors	Intel i5 or better	Intel i7 or better
Video Cards	Dedicated Graphics Card with 2 GB+ Memory NVidia Chipset re	Dedicated Graphics Card with 4 GB+ Memory commended
Memory (RAM)	Workstation: 8 GB or better Server: 16 GB or better	Workstation: 16 GB or better Server: 32 GB or better
Hard Drives	Workstation: 80+ GB Server: 500+ GB for data storage	Workstation: 80+ GB SSD Server: 1 TB for data storage
Operating Systems	Workstations: Windows 10, Windows 11 – Professional versions Check compatibility of hardware dr Servers: Windows Server versions: 2016 / 2019 / 2022 Windows OS versions should be M	ivers with version of Windows. Windows Server versions: 2016 / 2019 / 2022
Monitors	Widescreen monitor (1366 x 768 minimum resolution)	Widescreen monitor 1920 x 1080 resolution
Networking	1 Gigabit WIRED Network Card Wireless networking not recommended	
Image Capture Device		

OFFICE COMPUTER SPECIFICATIONS

Recommended

Processors	Intel i5 or better	
Operating Systems	Windows 10, Windows 11 – Professional versions	
Memory (RAM)	8 GB or better	
Hard Drives	80+ GB	
Monitors	s Widescreen monitor (1366 x 768 minimum resolution)	
Backup System	Backup System Windows approved backup system for CADI Patients' images folder Backups should be stored off site.	